



## **Synapse Mobility V7 IHE Integration Statement**

SYNAPSE is a registered trademark of FUJIFILM Medical Systems U.S.A., Inc. in the United States, Australia and Canada.

SYNAPSE is a registered trademark of FUJIFILM Corporation in China, Japan, Austria, Belgium, Denmark, Finland, France, Germany, Greece, Ireland, Italy, Luxembourg, The Netherlands, Portugal, Spain, Sweden, United Kingdom, Latvia, Slovakia, Slovenia, Cyprus, Lithuania, Hungary, Czech Republic, Malta, Poland, Estonia, Bulgaria, and Romania.

FUJIFILM Medical Systems U.S.A., Inc.  
81 Hartwell Avenue Suite 300  
Lexington, MA 02421  
Telephone: 1-800-272-8465

<b>IHE Integration Statement</b>		<b>Date:</b> July 5, 2019
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>
FUJIFILM Medical Systems U.S.A., Inc.	Synapse Mobility	V7.1
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>
Consistent Time (CT)	Time Client	None
Scheduled Workflow	Image Display	None
Consistent Presentation of Images (CPI)	Image Display	None
Key Image Note (KIN)	Image Display	None
Simple Image and Numeric Report (SINR)	Report Reader	None
Access to Radiology Information (ARI)	Image Display	None
Cross-Enterprise Document Sharing (XDS.b)	Document Consumer	None
Patient Demographics Query (PDQ)	Patient Demographics Consumer	None
<b>Internet address for vendor's IHE information:</b> <a href="http://www.fujifilmusa.com/products/medical/technical-documents/product-bulletins/">http://www.fujifilmusa.com/products/medical/technical-documents/product-bulletins/</a>		
<b>Links to Standards Conformance Statements for the Implementation</b>		
<b>DICOM</b>	<a href="http://www.fujifilmusa.com/products/medical/technical-documents/product-bulletins/">http://www.fujifilmusa.com/products/medical/technical-documents/product-bulletins/</a>	
<b>DICOMweb</b>	<a href="http://www.fujifilmusa.com/products/medical/technical-documents/product-bulletins/">http://www.fujifilmusa.com/products/medical/technical-documents/product-bulletins/</a>	
<b>Links to general information on IHE</b>		
North America: <a href="http://ihe.net/IHE_Worldwide/">http://ihe.net/IHE_Worldwide/</a>	Europe: <a href="http://ihe.net/IHE_Worldwide/">http://ihe.net/IHE_Worldwide/</a>	Asia/Pacific: <a href="http://ihe.net/IHE_Worldwide/">http://ihe.net/IHE_Worldwide/</a>